

Erlab 632 Ductless Fume Hood

Captair Smart medium ductless filtering fume hood



Features

- Erlab’s unique Flex Filtration can manage a large majority of laboratory applications
- Erlab Captair Smart fume hood are featured with configurable filtering systems above the enclosure, which can be adapted to filter gases, solvents, powders and particulates, keeping the user and lab protected and efficiently recirculating purified air back into the room, free of any toxic hazardous vapours or odours
- The smart technology of a simple, intuitive communication by pulsing light indicates the status of the fume hood, alerting the user when; air flow is compromised, there is a fan failure or filter breakthrough
- Intelligent side sensors detect when the filter is saturated and requires replacement
- Erlab’s smart technology offers low power consumption and allows the unit to operate with a standard plug using only 65W
- The 6° slanted sash offers optimal visibility and an ergonomic workspace all while the unit is automatically monitoring face velocity
- Large access opening on the front sash allows easy access for your work
- Full remote control and monitoring via the Erlab eGuard app allows users to change settings and delivers safety alerts immediately to a mobile, tablet or PC
- Erlab Captair Smart fume hoods meet full compliance with the international safety standards

Product Code	Product Description	Pack Size	Quote
FCC-632-SMART	Captair Smart 632 Ductless Fume Hood	ea	
FCA-M632-RWT(T)	Black Resin Worktop to suit 632	ea	
FCA-M632-GWT	Glass Worktop to suit 632	ea	
FCA-M632-CLEAR	Clear Transparent Back Panel to suit 632	ea	
FCA-M632-BENCAP	Benchcap Fixed Table with Utility Cannel to suit 632	ea	
FCA-M632-SHELF	Optional Shelf for Benchcap to suit 632	ea	

Technical Data

Dimensions (mm)	
Product Outer (H x W x D):	1106/1292 x 1600 x 615 mm
Air Flow (m3/h)	660
Filtration Technology	3
Internal Lighting	650
Power Consumption W	160

Relevant Product

- [321S Fume Hood](#)
- [321 Fume Hood](#)
- [391 Fume Hood](#)
- [392 Fume Hood](#)
- [481 Fume Hood](#)
- [483 Fume Hood](#)
- [633 Fume Hood](#)
- [714 Fume Hood](#)
- [Bio PCR Workstation](#)
- [Laminar Flow Fume Cabinet](#)
- [Pyramid Glove Box](#)
- [Safety Weighing Cabinets](#)
- [Green Fume Hood](#)
- [Customised Instrument Enclosures](#)
- [Ministore Vented Cabinet](#)
- [834 Filtration Storage Cabinets](#)
- [Filtration Boxes for Safety Cabinets](#)
- [Halo Air Purifier](#)

Filter Replacement

Do your filters in your Erelab Fume Hood need changing? Please click the button below, complete and submit the form and an Ecolab Life Science specialist will be in touch with a filter replacement quote.

- [Filter Replacement](#)

Applications

- Life Science
- Healthcare Facilities
- Hospitals

Resources

- [Video](#)
- [Video](#)
- [Brochure](#)